

This Patent Marking page is intended to provide information on U.S. Patents pursuant to 35 U.S.C. §287(a). The information on this page is effective from **May 1, 2025** until it is next revised. The column on the left side of the page identifies the product and the column on the right side of the page indicates the one or more patents that may cover that product.

PRODUCT	PATENT(S)
(H1 System) Tork Matic® Hand Towel Roll Dispensers	D648,585; D641,578 ; 8,561,441; 8,609,216; 8,733,571; 7,837,077; 8,955,790; 8,730,643; 8,490,913; 9,326,647; 9,918,597; 9,986,874; 10,099,412; 10,568,470 10,238,245; 10,987,839; 11,363,919 D811,120; 9,326,647; 10,136,768; 10,143,340; 10,674,876;
(H1 System) Tork® Refills	D716,149; D757,539; 8,006,929 7,975,957; 8,579,184; 10,093,467
(H2 System) Tork Xpress® Multifold Hand Towel Dispenser, (H23 Systems) Tork Xpress® Countertop Multifold Hand Towel Dispenser	D903,365; D903,366, 8,561,441; 8,609,216; 8,733,571; 10,099,412; 10,398,265; 10,674,876; 10,987,839; 10,136,768
(H2 System) Tork® Refills	D716,149; D757,539; 9,556,561; 8,579,184; 10,093,467; 11,753,225
(H5 System) Tork PeakServe® Continuous Hand Towel Dispenser; Elevation, Recessed	8,561,441; 8,609,216; 8,733,571; 9,999,325; 10,099,412; 10,058,221; 10,188,244; 10,568,471; 10,390,664; 10,758,096; 10,842,329; 10,987,839; 11,166,604; 11,206,956; 11,484,162; 11,812,899; 11,878,450; 11,871,877; 11,957,279; 11,963,642
(H5 System) Tork PeakServe® Refills	D716,149; D757,539; 8,877,313; 9,556,561; 9,004,313; 9,668,621; 10,093,467; 10,413,136; 10,875,705; 11,753,225

(H25 System) Tork® C-Fold hand towels	D716,149; D757,539; 8,579,184; 10,093,467
(H71 System) Tork® Hand Towel Roll Dispenser	11,179,011
(H71 System) Tork® Refills	D716,149; D757,539; 8,579,184; 10,093,467
(H76 System) Tork® Refills	8,579,184; 10,093,467
(H80 System) Tork® Mechanical Hand Towel Roll Dispenser	11,179,011
(H80 System) Tork® Refills	D716,149; D757,539; 8,579,184; 10,093,467
(H86 System) Tork Refills	D716,149; D757,539; 8,579,184; 10,093,467
(M1 System) Tork® Elevation Mini Centerfeed Hand Towel Dispenser, (M2 System) Tork® Elevation Centerfeed Hand Towel Dispenser	8,561,441; 8,609,216; 8,733,571; 9,456,717; 10,099,412; 10,987,839
(M1 System) Tork® Advanced Soft Mini Centerfeed Hand Towel Refill, (M2 System) Tork® Wiper Centerfeed Dispenser Refills	D716,149; D757,539; 8,318,281; 8,579,184; 10,093,467
(T1 System) Tork® Jumbo Bath Tissue Roll Dispenser; Elevation line	8,561,441; 8,609,216; 8,733,571; 10,099,412; 10,987,839
(T1 System) Tork® Refills	D716,149; D757,539; 8,579,184; 10,093,467
(T2 System) Tork® Image Design Mini Jumbo Bath Tissue Roll Dispenser; Elevation line, (T2 System) Tork® Mini Jumbo Bath Tissue Roll Dispenser; Elevation line, (T2 System) Tork® Twin Mini Jumbo Bath Tissue Roll Dispenser; Elevation line	D674,634; 8,561,441; 8,609,216; 8,733,571; 10,099,412; 10,136,768; 10,674,876; 10,987,839;
(T5/T7 System) Tork® Coreless High Capacity Toilet Paper Dispenser; Elevation line, Tork® Optiserve Coreless High Capacity Toilet Paper Dispenser; Elevation line	D788,486; D788,488; D788,487; D789,111; D674,634; 8,561,441; 8,609,216; 8,733,571; 7,513,453; 10,099,412; 10,987,839; 11,559,174; 11,897,715
(T7 System) Tork® Coreless High Capacity Spindle Kit Types A, B, C, D	D788,486; D788,488; D789,111, D788,487
(T10 System) Tork® High Capacity Bath Tissue Roll Dispenser for Opticore®	9,648,995;
(T10 System) Tork® High Capacity Bath Tissue Roll Dispenser for Opticore® Refills	8,579,184; 10,093,467

(T11 System) Tork® Toilet Paper Roll Dispenser for Opticore® (horizontal)	9,648,995
(T11 System) Tork® Refills	8,579,184; 10,093,467
(T24 System) Tork Bath Tissue Roll	D716,149; D757,539; 8,579,184; 10,093,467
(T26 System) Tork® High Capacity Bath Tissue Roll Dispenser; Elevation line	D788,486; D788,487; 7,967,235; 8,511,600; 8,561,441; 8,609,216; 8,733,571; 10,099,412; 10,987,839, 11,529,029
(T26 System) Tork® High Capacity Bath Tissue Roll Dispenser Refills	D716,149; 8,511,600; 8,579,184; 10,093,467
(T34 System) Tork® Premium Bath Tissue Roll	8,579,184; 10,093,467; D716,149; D757,539
(S1 System) Tork® Liquid Skincare Dispenser; Elevation line	D634,142; D662,743; 8,561,441; 8,609,216; 8,733,571; 8,887,964; 9,180,474; 8,701,942; 8,695,856; 8,708,200; 10,099,412; 10,987,839
(S1 System) Tork® Liquid Skincare Dispenser Refills	8,701,942; 8,695,856; 8,708,200; 10,470,616
(S4 System) Tork® Image Design Foam Skincare Manual Dispenser (S4 System) Tork® Image Design Foam Skincare Automatic Dispenser	D809,819; 8,561,441; 8,609,216; 8,733,571; 8,807,398; 10,136,768; 10,674,876
(S4 System) Tork® Elevation Foam Skincare Manual Dispenser (S4 System) Tork® Elevation Foam Skincare Automatic Dispenser with Intuition Sensor	D634,142; D662,743; 8,561,441; 8,609,216; 8,733,571; 8,807,398; 10,099,412; 10,987,839
(S4 System) Tork® Refills	D1,043,166; 8,807,398; 10,213,062; 11,051,660; 12,121,916
GFS – ARRAY® Extra Mild Foam Handwash refill	8,807,398; 10,213,062
(N4 System) Tork® Image Xpressnap® Napkin Dispenser, Aluminum , (N10 System) Tork® Image Xpressnap Café® Napkin Dispenser, Aluminum	D745,336; 9,986,876; 10,582,814
(N4 System) Tork® Image Xpressnap® Napkin Dispenser, Walnut, (N10 System) Tork® Image Xpressnap Café® Napkin Dispenser, Walnut	D745,337

(N4 System) Tork® Signature Xpressnap® Tabletop Napkin Dispenser, (N10 System) Tork® Classic Xpressnap Café® Napkin Dispenser	D705,016; D705,015; 7,178,689; 9,622,627; 9,901,222
(N4 System) Tork® Signature Xpressnap® Stand Napkin Dispenser, (N4 System) Tork® Classic Xpressnap® Stand Napkin Dispenser with Drive Thru Face Plate	D675,851; D701,434; 10,687,995, 8,328,046; 10,405,705
(N4 System) Tork® Signature Xpressnap® Counter Napkin Dispenser	D692,730; 9,622,627; 9,901,222; 9,986,876
(N4 System) Tork® Signature Xpressnap® In-Counter 5", 12", 20" Napkin Dispenser	D705,017; 9,433,331; 8,991,646; 8,863,984; 9,986,876; 10,582,814
(N4 System) Tork® Signature Xpressnap® Drive-Thru Dispenser	D740,049; D729,541; 10,405,705; 10,729,289; 10,687,995
(N4 System) Tork® Refills (napkins), (N10 System) Tork® Refills (napkins)	D716,149; D757,539; 8,579,184; 10,093,467
(N14 System) Tork® Xpressnap Fit Tabletop Napkin Dispenser , (N14 System) Tork® Xpressnap Fit Counter Napkin Dispenser –	D692,730; D705,015; D705,016; D995,123; 9,986,876; 9,622,627 9,986,876; 12,042,099; 12,144,470
(N14 System) Tork® Xpressnap Fit Dispenser Napkin Refills –	D915,080
Tork Paper Wipers (Refill ID: M2, W1, W2, W4, W24)	D719,842; D719,836; D719,460; D719,841; 8,318,281
Tork exelCLEAN (Refill ID: W1, W2, W4, W24, W7, W8)	D770,187; 8,906,816
(W2 System) Tork® Maxi Wiper Centerfeed Dispenser	8,561,441; 9,456,717
(W4 System) Tork® Folded Wiper Cloth Dispenser	8,561,441
(W6 System) Tork® Performance Washstation Dispenser	D775,849; D791,498; 8,561,441
(W6 System) Tork® Performance Washstation Dispenser Refills	8,318,281
(W8 System) Tork® Antimicrobial Foodservice Cloth	D720,234; wa1347,432,219; 8,318,281; 8,906,816
Tork® Wipers ¼ fold	D719,842; D719,836; D719,460; D719,841; 8,906,816; 8,318,281
Tork® Wipers Handy box	D719,842; D719,836; D719,460; D719,841
Tork Visison cleaning (former Easy Cube)	D732,530; 9,892,372; 10,390,365; 10,505,794; 10,973,377; 11,769,109; 11,849,886; 12,224,976

Tork LinStyle® Valet Guest Towel 1/6 Fold	9,303,362
Shamrock 68662 Dispenser	11,179,011
Swish 772CG	11,179,011
Renown RENOS5174WB and RENO5173WB	11,179,011
Buckeye 177286	11,179,011
Brady PH8682	11,179,011
Prolink 77279 and 77289	11,179,011
Boardwalk BWK29GREEN	11,179,011
Prime Source 75000882	11,179,011
Hillyard 432190 and 432690	11,179,011